



Active[®]

by ACCESS[™]



Removable shopping basket located under the seat. Basket contents are protected from theft and weather.



Handle height is easily adjustable with no loose hardware. Positions are clearly numbered for reference.



Four-position seat height (51-64cm) adjusts in seconds with dual release levers. No tools required.



Patented maintenance-free brake system features internal brake cables for user safety.



Curb climber makes it easy to get the Active over curbs and thresholds, even with a full shopping basket.



Patented indoor/outdoor ski solution makes the Active a year-round companion. Skis snap on in seconds with no screws or tools. Optional.



Effortless one-handed folding. Active stands and locks in the folded position.

Adapts with ease.

This type of rollator has never been seen before. Active has an attractive design, and stands out from the other rollators by combining all of the best user-friendly features into one and the same rollator. Both the handle and seat height can easily be adjusted with a simple pull of a lever, no tools required. The Active easily folds together for parking or transporting, and the shopping bag is conveniently located under the seat. Large diameter wheels make the Active a practical walking aid both indoors and out.

This rollator is the perfect companion for taking strolls outdoors, shopping trips to the local supermarket, or enhancing your safety in and around the house. You can't find a more user friendly or versatile walker anywhere in the world. You don't have to sit inside your home and watch the world go by.

Get Active, and stay active!

TECHNICAL SPECIFICATIONS

Type	Access Active
User height	160-190 cm
Max user weight	125 kg (275 lbs)
Dimensions (cm)	60 x 84 x 79-100
Folded dim. (cm)	60 x 50 x 89
Weight	8,9 kg (20 lbs)
Materials	Aluminium and Plastic
Cleaning	Detergent, steam or high pressure water.

Featuring a full 3 year manufacturers guarantee.

International Pending Patents: NO 20034379; US 10/895,903; CN 200410060031.2; PCT/N004/00275; NO 20040015; EP 000196134-0001; US 29/207,707; CN 200430046440.8

Access AS also designs, develops, and produces other products that can make your day easier. Contact us for more information.